




AGENTIC AI ARCHITECTURE DIAGRAMS


SYSTEM OVERVIEW (ASCII)










DATA FLOW DIAGRAM





STEP 2: Autonomous Analysis (Every 30s)





STEP 3: UI Updates (Real-time)

WebSocket Events:


'sensor_update' —————▶ Update patient cards
Update vital stats

'ai_alert' —————▶ Show notification
Play sound (critical)
Update alerts panel


'analysis_result' —————▶ Display AI summary
Show risk scores
Update agent status

MULTI-AGENT COLLABORATION





EXAMPLE SCENARIO: Fall Detection



+ BATHROOM 2 in 7 days" risk: 0.85"

CRITICAL Pattern Urgent



COORDINATOR:

"Patient fell in
bathroom. 2nd fall
this week. High
injury risk."



ALERT AGENT:

 CRITICAL ALERT

- Dispatch emergency
- Check for head injury
- Call family



Broadcast to ALL devices

PATIENT STATE LIFECYCLE

PATIENT STATE MANAGEMENT

NEW DEVICE DETECTED





```
Create PatientState |
{                  |
  device_id        |
  history: []      |
  alerts: []       |
  risk_score: 0    |
  patterns: []     |
}
```



CONTINUOUS DATA STREAM




```
┌───────────▼───────────┐
| Add data | ← MongoDB Listener
| to      | (Every 1s)
| history |
└──────────┘
```




```
History Buffer |
(Max 100 points) |
[              |
  {HR:75, SpO2:96}, |
  {HR:78, SpO2:95}, |
  ...           |
]
```





▼

CLEANUP (Optional)



STATE QUERY (On-demand)



Return current state:

```
{  
  "device_id": "ABC123",  
  "risk_score": 0.42,  
  "data_points": 1523,  
}
```

```
"last_update": "2025-01-30T15:23:45",  
"patterns": [...]  
}
```

ALERT WORKFLOW

ALERT GENERATION WORKFLOW

TRIGGER: Anomaly Detected



```
Monitor Agent |  
returns: |  
{ |  
  severity: |  
  "CRITICAL", |  
  anomalies: [ |  
    "FALL", |  
    "BATHROOM" |  
  ] |  
}
```



```
Alert Agent |  
- Classify priority |
```

- Generate message

- Suggest actions



Alert Object


```
{
  id: "ALERT_12345",
  timestamp: "...",
  severity: "CRITICAL",
  message: "🔥 FALL...",
  actions_required: [
    "Dispatch EMT",
    "Check injury"
  ],
  auto_notify: true
}
```



Save WebSocket External

to Emit Notify

alerts[] → (SMS/Email)



UI Notification

- Show red banner

- Play sound

- Log to alerts panel
- Update agent status

🔧 SYSTEM STARTUP SEQUENCE

SYSTEM INITIALIZATION

START: python agentic_medicore.py



1. Import modules

- Flask
- SocketIO
- pymongo
- aisuite



2. Initialize

- AI_CLIENT
- PATIENT_STATES
- ACTIVE_ALERTS



3. Start threads

Thread 1:

MongoDB

Listener

Thread 2:

Autonomous

Monitor

4. Start Flask

- Load routes

- Bind SocketIO

- Listen 5000

5. Ready to serve

✓ Agents ON

✓ MongoDB ON

✓ WebSocket ON

RUNTIME LOOPS (Concurrent):

Loop 1: MongoDB Listener

while True:

 Check new documents (1s interval)

 Update patient states

 Emit sensor_update

Loop 2: Autonomous Monitor

while True:

 For each patient:

 Run Monitor Agent

 Run Analyzer Agent

 Run Predictor Agent

 Coordinate results

 Generate alerts if needed

 Sleep 30s

Loop 3: WebSocket Server

Event-driven:

 On connect → Send current states

 On request_analysis → Run analysis

 On disconnect → Cleanup

CLIENT-SERVER COMMUNICATION


WEBSOCKET COMMUNICATION FLOW

CLIENT SIDE (Browser)

|
| socket.io.js





REAL-TIME UPDATES:





DEPLOYMENT ARCHITECTURE





Gunakan diagram ini untuk memahami alur kerja system! 🚀